

# PROMICRA

## QuickPHOTO MICRO 3.1

System for Digital Photomicrography, Image Editing and Measurements Designed for Biological and Life Science Applications

The screenshot displays the QuickPHOTO MICRO 3.1 software interface. The main window shows a grid of images with measurement overlays. The top-left image shows a wasp with a green line segment labeled  $L1=13.64\text{ mm}$  and a 2 mm scale bar. The top-right image is a close-up of a wasp's wing. The bottom-left image shows a textured surface with a 0.2 mm scale bar. The bottom-right image shows a circular object with a red outline, labeled  $P4=4.47\text{ mm}$  and  $A4=1.54\text{ mm}^2$ . A 'Table of Measured Values' window is open at the bottom, displaying the following data:

Object Type	Measurand	Meas. Nr.	Value	Units	Description	Measurand	Length
Line Segment	Length	1	13.65	mm		Number of Measurements	1
						Minimal Value	13.65
						Maximal Value	13.65
						Average Value	13.650
						Standard Deviation	0.000

The interface includes a menu bar (File, Edit, Image, View, Measurements, Annotations, Imaging Devices, Modules, Reports, Window, Help), a toolbar, a Thumbnail Browser on the left, and a status bar at the bottom showing '2550 x 1919 (19.98mm x 15.04mm) Lumenera® 8 INFINITY Digital Camera'.

# QuickPHOTO MICRO 3.1

QuickPHOTO MICRO 3.1 program is designed for acquiring of digital images from microscopes equipped with digital cameras, editing, annotating and saving of acquired images and measurements. This program is designed for use especially in biological and life science application.


QuickPHOTO MICRO 3.1 program is equipped with many measuring functions: measuring of lengths, perimeters, areas, angles, counting of objects, phase analysis. Using the function for time-lapse image acquisition, it is possible to acquire images in a predefined time interval.

## Main Features Overview:

- **Image acquisition from a wide range of imaging devices** (microscopic digital cameras, DSLR cameras, analog cameras, scanners)
  - **Live view on a computer screen**
  - **Measurements in acquired images:** insertion of a calibrated scale bar, measuring of dimensions (lengths), perimeters, areas, angles, objects counting, phase analysis, overlay orthogonal grid display
  - **Measurements in a live view\*:** calibrated scale bar, dimensions (lengths), perimeters, areas, angles, counting of objects, overlay orthogonal grid display
  - **Measured values are displayed in a well-arranged table with basic statistics.** Particular measurements can be annotated. Measured values can be exported (optionally with the image with measuring objects) to a Microsoft® Excel® file or a text file.
  - **Advanced functions for saving of multiple images at once**
  - **Option to automatically save acquired images to a predefined folder**
  - **Microsoft® Word reports**
  - Support for zooming units, magnification changer units and magnification encoders (devices for automatic determination of currently used microscope magnification)
  - **Image editing functions, indication and annotation of interesting areas of the image, image descriptions**
  - **Advanced printing functions** including interactive preview
  - **Automatic time-lapse image capturing** with the options to control computer controlled illuminators and time-lapse HD video creation
  - Option to acquire images using a remote trigger (optional accessory)
  - Calibrations can be locked with a password protection; the calibrations are protected by a digital signature against unauthorized modifications
  - Support for work on two-monitor workstations
  - **Multilingual user friendly graphical interface**
  - Possibility to extend functions with additional modules (Deep Focus, RECORD IT, FLUO+)
- \*These features are available only for some of the supported imaging devices.

## Measuring Functions (in acquired images):


- **Insertion of a calibrated scale bar** (options to define type, size and colors; vertical or horizontal position)
- **Measuring of dimensions** (line segment and polyline tools)
- **Measuring of areas and perimeters** (ellipse, rectangle and polygon tools)
- **Measuring of angles** (arbitrary angle, angle from the X axis, angle from the Y axis)
- **Counting of objects**
- **Phase analysis**
- **Display of a calibrated overlay orthogonal grid**
- Options to define measurands and their labeling for individual measuring tools
- Measured values can be exported (optionally with the image with measuring objects) to a Microsoft® Excel® file or a text file.



# QuickPHOTO MICRO 3.1


## Measurements in a Live View\*:

- Measurements directly in a live view (no need to acquire an image to perform measurements)
  - **Insertion of a calibrated scale bar** (options to define type, size and colors; vertical or horizontal position)
  - **Measuring of dimensions** (line segment and polyline tools)
  - **Measuring of areas and perimeters** (ellipse, rectangle and polygon tools)
  - **Measuring of angles** (arbitrary angle, angle from the X axis, angle from the Y axis)
  - **Counting of objects**
  - **Display of a calibrated overlay orthogonal grid**
  - Measured values can be exported to a Microsoft® Excel® file or a text file
- \*These features are available only for some of the supported imaging devices.



## Phase Analysis

This feature enables measuring of percentage amount of a particular phase (color/shade of grey) in the image. Up to four phases can be measured in one image. It is possible to select various colors of displayed phases, add another color/shade of grey to the current phase, invert the phase, set tolerance and overlays. As a result the percentage amount of each phase in the image is calculated.



## Time-Lapse Image Capturing

The built-in *Timer* module enables a time-lapse image capturing (automatic image capturing in defined time interval). The *Timer* module is capable to control the computer controlled illuminators. That way the illumination can always be switched on only for the time necessary for image acquisition to prevent a negative long-term effect of the illumination to the specimen.

## Time-Lapse Video Creation


The built-in *Time-Lapse Video* module can be used to create HD video clips (with resolution up to 1920 x 1080 pixels) from time-lapse captured images. For example when several (or several thousands) images of a slowly moving or changing object are acquired, the *Time-Lapse Video* module can be used to create a video file in which the movement will be recorded in a speeded-up form.

## Image Editing Functions

The images are immediately after acquisition transferred to a computer and prepared for processing. QuickPHOTO MICRO 3.1 program enables work with multiple images at the same time. For comparison two or four images can be displayed side by side. The *Thumbnail Browser* enables easy switching to another opened image by a single mouse click. The *Navigator* window facilitates the work with the magnified images.

The images can be edited using the following functions:

- **Crop** (crop size can be defined in pixels,  $\mu\text{m}$  or mm)
- **Rotation** (by 90° clockwise and counter clockwise)
- **Flip and mirror**
- **Resize**
- **Brightness and contrast adjustment**
- **Color balance adjustment**
- **Histogram adjustment**
- **Filters** (sharpen, blur, negative, grey scale, median filter)
- **White balance and black balance**
- **Indication of interesting areas** (line, arrow, ellipse, rectangle, circle and freehand drawing tools)
- **Insertion of text annotations**



# QuickPHOTO MICRO 3.1

## The Advantage of Working with Layers

All inserted objects (drawings, annotations, measurements) are displayed in a separate layer “above the image“. It is possible to modify (change of the line's color and thickness, change of size, change of position), delete, or hide them at any time. It is also possible to delete the whole layer with inserted objects or merge it with the image.

## Microsoft® Word Reports


Images with measured values can be easily exported into reports in Microsoft® Word format. The user can create custom report templates with various layouts and number of images. Templates can be created in Microsoft® Word program using standard formatting functions.

## Control Panels of Imaging Devices

Each of the supported imaging devices can be controlled by a dedicated control panel. The control panel contains a live view from the imaging device and elements for control of its functions and properties.

## Imaging Devices Controlled via TWAIN Interface

The imaging devices, which do not have their own control panel in QuickPHOTO MICRO 3.1 program, can be controlled via TWAIN interface. Those devices can also be used in conjunction with built-in *Timer* module as well as with *Deep Focus* z-stacking additional module.



## Image Acquisition from Analog Video Sources

The QuickPHOTO MICRO 3.1 program is equipped also with a control panel for image acquisition from analog video sources (analog cameras, DVD players and other analog devices).

## Two-Monitor Workstations Support

Two-monitor workstations are supported by QuickPHOTO MICRO 3.0 program. The main program window can be displayed on a primary monitor, while the control panel of a digital imaging device with a live view from a microscope can be displayed on a secondary monitor.

## Remote Trigger

A remote trigger is an optional accessory which makes image acquisition even more comfortable. The images can be acquired by pressing a trigger button.

## Minimal System Requirements:

Processor	Operating Memory	USB Ports	Operating System	Minimal Monitor Resolution
Single-core 2.4 GHz or multi-core	1 GB	2x USB 2.0/3.0	Microsoft® Windows® XP(SP3)/Vista/7/8/10	1024 x 768 pixels

## Supported Imaging Devices:

Complete list of currently supported devices can be found on the website: [www.promicra.com](http://www.promicra.com).

[www.promicra.com](http://www.promicra.com)