# **PRO-MSILL (SUNFLOWER)**

**[EN] User Guide**

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# **PRO-MSILL (SUNFLOWER)**

Copyright © PROMICRA, s.r.o. 2019

---

# Table of Contents

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
---

# [EN] PRO-MSILL (SUNFLOWER) - Universal 6-arm Diffuse LED Illuminator for Stereomicroscopes

## 1. Preface

PRO-MSILL (SUNFLOWER) is an universal 6-arm diffuse high power LED illuminator for stereomicroscopes. The illuminator consists of two parts:

- Illuminator unit PRO-MSILL
- Control unit PRO-MSILL-CB

## 2. Safety Precautions

### General

1. Be sure to set the main switch to “O” (OFF) before connecting or disconnecting the connection cables and power cord. To prevent electric shock and malfunction, never disconnect or connect a connection cable while the main switch is in “I” (ON) position.
2. Always use the specified AC adapter only. Using another AC adapter may cause the device to malfunction as well as catch fire or burn due to abnormal heating.
3. The cords and cables are vulnerable to bending or twisting. Do not apply excessive force to them and do not bend them severely. Lay out the cables so that they are not in contact with the heat generating devices.

### Computer and Software

The data on the computer may be destroyed by an unexpected event. Be sure to keep a backup of the data.

1. PROMICRA will not assume any liabilities for any damage incurred due to the use or no-usability of this system, including compensation for the loss of data.
2. The computer used with this system should be set up with and run Microsoft® Windows® 10/8.1/8/7/Vista/XP.
3. For the OS in the computer, the user is requested to create a backup and retain it. For details on the computer and Microsoft® Windows® OS please refer to their respective manuals.
4. PROMICRA guarantees the quality of this product in factory shipment conditions. PROMICRA will not assume any liabilities for operation errors and functional faults incurred due to the alteration of the environmental setup (including BIOS changes), installation of other software or the addition of hardware to the computer by the user.
5. When the HDD free space reduces, the data processing speed may slow down or errors may occur frequently. To prevent this, delete unnecessary data files frequently. For how to delete data files, refer to the manual for Microsoft® Windows® OS.

6. Do not open the enclosure of the computer and touch the power supply or the circuit board's heat generating sections right after use as it may burn your hand. Wait until the internal temperature drops sufficiently.
7. Sharp edges inside the computer may cut your fingers, so take extra care.

Use a computer that complies with safety standards of your country.

“CE” mark indicates that this product complies with corresponding European EN standards.

This symbol (crossed-out wheeled bin WEEE Annex IV) indicates separate collection of waste electrical and electronic equipment in the EU countries. Please do not throw the equipment into domestic refuse. Please use the return and collection systems available in your country for disposal of this product.


*RoHS compliance* – this product complies with the European Union's RoHS Directive, which restricts the use of the hazardous substances in electrical and electronic equipment.

### **3. Device Description - PRO-MSILL Illuminator**

The PRO-MSILL illuminator has the high power LEDs mounted on six tilting arms and a central ring. Each tilting arm contains two independently controllable LED segments. Optimal illumination conditions can be achieved by arbitrary combination of switched ON/OFF LED segments as well as by tilting the arms. 12+1 LED segments can be controlled. There are plastic diffusers mounted in front of each LED segment. The illuminator is equipped with mounting mechanism, which allows mounting to the objectives with diameter up to 69 mm.

[EN] PRO-MSILL (SUNFLOWER) - Universal 6-arm Diffuse LED Illuminator for Stereomicroscopes


---

## 4. Device Description - PRO-MSILL-CB Control Unit

### Top Side


Top side of the PRO-MSILL-CB control unit contains following elements:

1. Buttons to switch ON/OFF each LED segment of outer ring
2. Buttons to switch ON/OFF each LED segment of inner ring
3. Button to switch ON/OFF central LED ring
4. Button for rotation of current combination of switched ON/OFF LED segments counter-clockwise
5. Button for rotation of current combination of switched ON/OFF LED segments clockwise
6. Button to switch ON/OFF the whole inner ring of LED segments

7. Button to switch ON/OFF the whole outer ring of LED segments
8. Button for increasing the intensity of all LED segments
9. Button for decreasing the intensity of all LED segments
10. Button to switch ON/OFF the whole illuminator (standby mode)

## Rear Side

Rear side contains the following elements:

- *PWR* - main switch
- *DC* - connector for AC adaptor
- *USB* - connector for PC connection
- *ILL.* - connector for PRO-MSILL illuminator

## 5. Device Installation

### 5.1. Cabling

#### **Important**

The device is supplied with 24V AC adaptor. Do not use other AC adaptors, otherwise the damage of the device may occur.

When connecting cables, follow these steps:

1. First connect the PRO-MSILL illuminator cable to ILL. connector.

2. Then connect the AC adaptor to DC connector.

**Important**

Never connect the PRO-MSILL illuminator cable when the PRO-MSILL-CB control unit is switched on.

## 5.2. Mounting of PRO-MSILL Illuminator to a Microscope

1. Raise the microscope using the focus knob or remove the microscope from the stand.
2. Rotate the handle of the PRO-MSILL mounting mechanism in direction of the OPEN arrow to the final position.
3. Put the illuminator on the stage and point the tilting arm with white dot to the front. The cable should be directed towards the rear.
4. Place the microscope back to the operating position. Carefully insert the objective through the hole in the upper part of the illuminator.

**Important**


Do not use excessive force when manipulating the mounting mechanism of the illuminator.

## 5.3. PC Control Installation

The device can be controlled also by a computer. Install the device drivers following these steps:

1. Browse the installation disk, locate the *Drivers* folder and run the *CDM vx.xx.xx WHQL Certified.exe* file and follow the on-screen instructions until finished.
2. Connect the PRO-MSILL-CB control unit to a running PC with Microsoft® Windows® 10/8.1/8/7/Vista/XP operating system using the supplied USB cable.
3. Turn the PRO-MSILL-CB control unit on using the PWR switch at the rear side.
4. Wait until the device drivers will be automatically installed.

When the device driver has been properly installed “USB Serial Port (COM#)” should be displayed in Device Manager as shown below.

After this procedure the PRO-MSILL illuminator is ready to communicate with the PC.

## 6. Device Operation

### 6.1. PRO-MSILL Illuminator Operation

#### Important

Be careful, the illuminator unit may be very hot!

The illuminator can be used in two general ways:

1. The PRO-MSILL illuminator can be placed on the microscope stage and the objective can move freely through the illuminator's hole when changing the focus of the microscope.
2. The PRO-MSILL illuminator can be tightened up by the mounting mechanism to the objective. The illuminator will move up and down with the objective when changing the focus of the microscope.

#### Note

The illuminator provides best illumination conditions when fully closed.

### Placing a Sample under the Microscope

Raise one or more tilting arms using the plastic handles and place a sample under the microscope. Close the illuminator again.

**Note**

Before raising more than two tilting arms, it is recommended to tighten up the illuminator to the objective using the mounting mechanism.

## **Oblique Illumination**

To set an oblique illumination, tighten up the illuminator to the objective and raise the desired number of tilting arms using the plastic handles to desired angles.

## **6.2. PRO-MSILL-CB Control Unit Operation**

### **Powering up the Illuminator**

Power up the illuminator by switching the PWR switch to I (ON) position. All LED segment will be switched on to half intensity.

### **Setting the Illumination Intensity**

The intensity of all switched on LED segments can be adjusted in 32 steps by repetitive press of buttons 8 and 9. The intensity will change continuously while button 8 or 9 is being pressed. Current illumination intensity is signalled by one of 4 blue LEDs placed next to the buttons 8 and 9. When maximal illumination intensity is reached, the top blue LED will start flashing and the beep will sound. When the minimal intensity is reached the bottom blue LED will start flashing.

### **Switching the LED Segments ON/OFF**

Each one of the 12+1 LED segments can be independently switched ON/OFF by the buttons 1, 2 and 3 to achieve desired illumination conditions under the microscope. Those buttons layout corresponds to the positions of the controlled LED segments on the PRO-MSILL illuminator. The front tilting arm of the PRO-MSILL illuminator as well as the corresponding buttons of the PRO-MSILL-CB control unit are marked by a dot. The LED segments which are currently switched on are signalled by blue LEDs inside the corresponding buttons.

## **Switching ON/OFF the Whole Inner/Outer Ring of LED Segments**

The LED segments mounted on the tilting arms create two rings, the inner one and the outer one. To switch ON/OFF the whole ring of LED segments use buttons 6 and 7.

## **Rotation of the Current Combination of Switched ON/OFF LED Segments**

Current combination of switched ON/OFF LED segments can be rotated counter-clockwise/clockwise using the buttons 4 and 5.

# **7. Specifications**

## **PRO-MSILL Illuminator**

Dimensions	195 x 195 x 133 mm
Weight	1950 g
Number of independent LED segments	12+1
Number of LEDs	432
LED's lifetime	> 20 000 h
Compatible diameters of objectives	up to 69 mm
Environmental conditions	In use: temperature: 0 - 40 °C, humidity: 30 % - 80 %; In transportation: temperature: -25 - 60 °C, humidity: 10 % - 90 %
Manufacturer	PROMICRA, s.r.o., Czech Republic, <a href="http://www.promicra.com">www.promicra.com</a>

## PRO-MSILL-CB Control Unit

Dimensions	154 x 188 x 60 mm
Weight	635 g
PC interface	USB
AC adaptor	DC 24 V  1.25 A
Connectors	DC power socket, USB (B-type) - PC interface, 20-pin connector for illuminator
Environmental conditions	In use: temperature: 0 - 40 °C, humidity: 30 % - 80 %; In transportation: temperature: -25 - 60 °C, humidity: 10 % - 90 %
Supported operating systems	Microsoft® Windows® 10/8.1/8/7/Vista/XP
Manufacturer	PROMICRA, s.r.o., Czech Republic, <a href="http://www.promicra.com">www.promicra.com</a>

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