

# CB-ZM

---

## Control Unit for Focus Motorization of Optical Microscopes



CB-ZM control unit with MFD focus drive mounted on a stereomicroscope

### Features Overview:

- Suitable for focus motorization of almost all types of optical microscopes and stereomicroscopes
- Can be used separately or in conjunction with Deep Focus\* z-stacking module
- Manual focus control using UP and DOWN buttons
- SPEED button for selection of focus speed
- Current focus speed signalized by LEDs on the front panel
- Remote control by a two-button foot switch (optional)
- Possibility of PC control via USB interface (supported by Deep Focus module)
- Dimensions: 10.5 cm x 8.5 cm x 3.5 cm.

\*Deep Focus is an additional software z-stacking module for QuickPHOTO programs.