


Image Stitching Module

Live Manual Image Stitching Module for QuickPHOTO Programs



Features Description:

- Seamless high quality image stitching using live view of a camera and a manual stage
- Can be used for stitching of just a few images as well as for stitching of whole samples with resolution up to several gigapixels.
- Excellent illumination homogeneity using shading (background) correction feature
- Enables precise measurements in stitched images.
- Can be used with any type of optical microscope with manual X-Y stage.
- Suitable for biological as well as industrial applications
- Can be combined with HDR module.
- Processing on a hard drive minimizes memory consumption.

Required Equipment :

- QuickPHOTO CAMERA, QuickPHOTO MICRO or QuickPHOTO INDUSTRIAL program in version 3.0 or 3.1 (recommended)
- PROMICAM USB 3.0 camera
- Microscope with manual X-Y stage
- Computer with Microsoft® Windows® 7/8.1/10 operating system
- SSD drive (not necessary, but recommended for speeding up work when stitching larger number of images)

PROMICRA, s.r.o. – Evropska 39 – 160 00 Prague 6 – Czech Republic

Phone: +420 222 952 667 – Email: promicra@promicra.com

www.promicra.com