

# PRO-MLS

## Universal Heavy-Duty Motorized Stand for Microscopic and Macroscopic Applications



PRO-MLS stand with PRO-ZBAD adapter and CB-ZM control unit

### Features Overview:

- Focusing and z-movement by UP and DOWN buttons, change of focus speed using SPEED button
- Compatibility with [Deep Focus](#)\* z-stacking software module
- Adapters for a wide range of optical devices (stereomicroscopes, macro cameras...) and accessories (microscopic stages, LED illuminators...)
- High bearing capacity (up to 16kg) and rigidity, silent running
- Travel range 19.5 cm
- Remote control by a two-button foot switch (optional)
- Z-movement range protected by built-in limit switches

(\*) Additional software module for [QuickPHOTO](#) programs



PRO-MLS stand with PRO-FUAD adapter and CB-ZM control unit

### Example of Use:



Motorization of OLYMPUS® SZX16 stereomicroscope

OLYMPUS® is a registered trademark of OLYMPUS Corp.

PROMICRA, s.r.o. – Evropska 39 – 160 00 Prague 6 – Czech Republic  
Phone: +420 222 952 667 – Fax: +420 222 968 963 – e-mail: [promicra@promicra.com](mailto:promicra@promicra.com)

[www.promicra.com](http://www.promicra.com)