

PRO-SZM1

Additional Module for Focus Motorization of Stereomicroscopes



PRO-SZM1 module on SZ2-ST stand with CB-ZM control unit

Functions Overview:

- Focusing and z-movement by UP and DOWN buttons, change of focus speed using a SPEED button.
- Manual focus function maintained.
- Compatibility with [Deep Focus](#)* z-stacking software module (module for creation of images with extreme depth of focus).
- Option to attach LED illuminators and other accessories.
- Remote control by a two-button foot switch (optional).
- Z-movement range protected by built-in limit switches.

Compatible Stands:

- OLYMPUS® SZ2-ST.

(*) Additional software module for [QuickPHOTO](#) software

Examples of Use:

Motorization of OLYMPUS® SZ61 stereomicroscope



Motorization of OLYMPUS® SZX7 stereomicroscope

OLYMPUS® is a registered trade mark of OLYMPUS Corp.