


# RECORD IT Module

## Software Module for QuickPHOTO Programs for Real-Time Video Recording from Digital and Analog Imaging Devices

RECORD IT is an additional software module for QuickPHOTO programs. It is intended for real-time video recording from digital and analog imaging devices.



### Functions Overview:

- Real-time video recording from several types of digital cameras.
- Real-time video recording from analog cameras and other sources of analog video signal (e.g. from video endoscopic systems, ultrasonographs, VCRs, DVDs).
- Two-button foot switch can be used for video recording and image acquisition control.
- Automatic file name creation in case of multiple recordings.
- Automatic saving of created video files into selected folder.

### List of supported imaging devices can be found on:

<http://www.promicra.com/products-recordit.php>

### System Requirements:

- QuickPHOTO CAMERA, MICRO, Industrial in version 2.2 or higher.

PROMICRA, s.r.o. – Evropska 39 – 160 00 Prague 6 – Czech Republic  
Phone: +420 222 952 667 – Fax: +420 222 968 963 – e-mail: [promicra@promicra.cz](mailto:promicra@promicra.cz)  
[www.promicra.com](http://www.promicra.com)