

CB-ZM

Control Unit for Focus Motorization of Optical Microscopes



CB-ZM control unit with MFD focus drive mounted on a stereomicroscope

Features Overview:

- Suitable for focus motorization of almost all types of optical microscopes and stereomicroscopes
- Focus control using UP and DOWN buttons
- Selection of focus speed – 4 predefined speed levels available
- Current focus speed signalized by LEDs on the front panel
- Automatic speed incrementation when holding UP or DOWN button
- Remote control by a two-button foot switch (optional)
- Compact design (control unit and a hand switch in a single unit) saves space on your desktop
- PC control option via USB interface (supported by Deep Focus* module)
- Dimensions: 10.5 cm x 8.5 cm x 3.5 cm.

*Deep Focus is an additional software z-stacking module for QuickPHOTO programs.