


FLUO+ Module

Software Module for Multi-Color Fluorescence Imaging

FLUO+ is an additional multi-color fluorescence imaging software module for QuickPHOTO programs. Up to 8 digital images acquired with different fluorescence filters can be composed into one multi-color image.



Functions Overview:

- Up to 8 digital images acquired with different fluorescence filters can be composed into one multi-color image
- Histogram and intensity of each of the fluorescence images can be adjusted independently
- X and Y shift correction for each image
- Color of each of the fluorescence images can be user defined
- Module settings can be saved to a file and loaded at any time

System Requirements:

- QuickPHOTO CAMERA, QuickPHOTO MICRO, QuickPHOTO INDUSTRIAL in version 2.2 or higher

