


# Image Stitching Module

## Image Stitching Module for QuickPHOTO Programs



### Features Description:

- Seamless high-quality live image stitching using live view of a camera and a manual stage
- Stitching of images acquired by any imaging device
- Can be used for stitching of just a few images as well as for stitching of whole samples with resolution up to several gigapixels.
- Excellent illumination homogeneity using shading (background) correction feature
- Enables precise measurements in stitched images.
- Can be used with any type of optical microscope with manual X-Y stage.
- Suitable for biological as well as industrial applications
- Can be combined with HDR module or Deep Focus module.
- Processing on a hard drive minimizes memory consumption.

### Required Equipment :

- QuickPHOTO CAMERA, QuickPHOTO MICRO or QuickPHOTO INDUSTRIAL program in version 3.0 or newer
- Microscope with manual X-Y stage
- Computer with Microsoft® Windows® 7/8.1/10 operating system
- SSD drive (not necessary, but recommended for speeding up work, when stitching larger number of images)

PROMICRA, s.r.o. – Evropska 39 – 160 00 Prague 6 – Czech Republic

Phone: +420 222 952 667 – Email: [promicra@promicra.com](mailto:promicra@promicra.com)

[www.promicra.com](http://www.promicra.com)