

# PRO-SZM1

---

## Motorized Focus Drive for Stereomicroscopes with 76mm Mount



PRO-SZM1 motorized focus drive on SZ2-ST stand with CB-ZM control unit

### Features Overview:

- Motorized Z axis (travel and focusing)
- Controlled by CB-ZM control unit: focusing and travel by UP and DOWN buttons, change of focus speed by SPEED button, PC control option
- Manual focus function maintained
- Minimal step 65  $\mu\text{m}$
- Travel range 12 cm
- Travel range protected by built-in limit switches
- Option to attach LED illuminators and other accessories
- Remote control by a two-button foot switch (optional)
- Compatible with Deep Focus extended depth of field (EDF) software module for QuickPHOTO software

### Compatible Stands:

- OLYMPUS® SZ2-ST

## Examples of Use:



Motorization of OLYMPUS® SZ61 stereomicroscope



Motorization of OLYMPUS® SZX7 stereomicroscope

OLYMPUS® is a registered trade mark of OLYMPUS Corp.