

PRO-LM-LED-FLUO

Fluorescence LED Light Sources



PRO-LM-LED-FLUO LED Light Sources



Compatible Control Units: PR-LC30 (left) and PC Controllable Unit PR-LC22 (right)

Modern Sources of Excitation Illumination which Offer Following Advantages Compared to Standard Lamps:

- Lifetime over 20 000 hours
- Low power consumption and heat generation
- Smooth intensity adjustment
- PC control option

Available Variants:

- Blue (FITC)
- Green (TRITC)
- UV 365 nm, UV 400nm (DAPI)

Compatibility with Existing Fluorescence Illuminators for Following Microscopes:

- OLYMPUS® BX/CX/CKX series

OLYMPUS® is a registered trade mark of OLYMPUS Corp.